STAT/S® Handpieces - Quick Reference Guide

Please consult Operator's Manual before using the devise(s).

Installation

Checking the pressure readings:

1. Place a handpiece coupling onto an air pressure gauge, if available.
2. Place the gauge onto the dental unit tubing.
3. Shut off the water and place your foot onto the rheostat.
4. Adjust the drive air, chip air & water on the unit, as necessary, to reach the recommended pressure readings.

Pressure readings for the STAT/S SL & ML:

- Drive air = 2.1-4.2 bar (30-60 p.s.i.)
- Recommendation = 2.8 bar (40 p.s.i.)
- Return air < 0.5 bar (7 p.s.i.)
- Water = 0.8-2.5 bar (11-36 p.s.i.)
- Spray air = 1.0-4.0 bar (14-57 p.s.i.)
 - The minimum amount of water required for spray cooling is 50 mL per min. An insufficient amount of water may cause overheating of the tooth.

Checking the STAT/S SL & ML coupling:

All 5 O-rings must be in place and should not be damaged. Missing or damaged O-rings should be replaced, or they may result in premature failure or negatively affect handpiece performance.

Important notes for burs & diamonds:

Diamond and carbide use: Only use burs that have shafts compliant with ISO 1797-1 type 3, made from steel or tungsten carbide and that meet the following criteria:

- Shaft diameter from 1.59 -1.60 mm.
- Overall length up to max. 25 mm.
- Shaft clamping length of at least 11 mm.
- Cutting edge diameter up to max. 2 mm.
- Short shank burs can be used with the STAT/S SL 101.1 model.

Did you select the right coupler for your tubing?

Connex LED S10010

6 pin, MULTIflex style coupling with LED bulb, anti retraction valve.

Connex L S10002

6 pin, coupling with light connection, anti retraction valve.

Connex GR S11002

5 hole, coupling with fiber optic glass rod, anti retraction valve.

Manufactured by:
SciCan GmbH
Wangener Strasse 78 Tel.: +49 (0)7561 98343 - 0
88299 Leutkirch GERMANY Fax: +49 (0)7561 98343 - 699

SciCan | Your Infection Control Specialist™

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