

OPTIM® Wipes Quick Reference Guide

 1-800-667-7733

The product's label should be consulted before using this surface disinfectant.

OPTIM Ready-to-Use Wipes

- 1 Pull out OPTIM towelette at a 45 degree angle. Close lid.



- 2 If surface is not heavily soiled, clean and disinfect surface with one wipe and allow surface to remain wet for the contact time indicated on the label.



- 3 If surface is visibly soiled, clean surface with one towelette and discard. Use a second towelette to disinfect and allow surface to remain wet for the contact time indicated on the label.



OPTIM SOLUTION Liquid Spray

- 1 Spray Ready-to-Use OPTIM liquid onto a SciCan Dry Towelette. Do not directly spray onto the surface.



- 2 If surface is not heavily soiled, clean and disinfect surface with one wipe and allow surface to remain wet for the contact time indicated on the label.



- 3 If surface is visibly soiled, clean surface with one towelette and discard. Use a second towelette to disinfect and allow surface to remain wet for the contact time indicated on the label.



IMPORTANT:
Do not allow OPTIM to pool on surfaces. Use of gauze dipped in solutions for wiping may cause pooling and should not be used. Do not spray directly onto the surface.



FOR BEST RESULTS:

- Do not over apply OPTIM.
- Use of OPTIM on brass, copper, anodized aluminum, carbon tipped instruments (eg. burs), and soft plastic; specifically, VITA SHADE, ROOT ZX APEX LOCATOR'S SCREEN and PELTON & CRANE CHAIR CONTROLS could discolour/or damage surfaces, if not used as directed and over applied.
- Use of gloves with corn starch may cause yellowing of surfaces.
- Use of gloves produced with high levels of sulphur may result in a foul odour.
- Dental Association guidelines recommend the use of utility gloves when cleaning and disinfecting surfaces.

OPTIM is safe, it does not require any GHS (Globally Harmonized System) hazard pictograms, signal words or statements.

PLEASE REFER TO THE ORIGINAL CONTAINER FOR COMPLETE DIRECTIONS AND ADDITIONAL INFORMATION.

Procedure for disinfecting PVS elastomeric impression materials with OPTIM

1. Once the impression is taken, rinse it well under running tap water to remove blood and saliva.
2. After rinsing, gently shake the impression within the sink basin to remove any adherent water with minimal splatter.
3. Spray the entire impression with OPTIM. Use the contact times as stated on the label.
4. Rinse the impression thoroughly with tap water to remove any excess disinfectant. After a thorough rinse, gently shake the impression inside the sink basin to remove any adherent water with minimal splatter.
5. Once dry, place the disinfected impression in a sealed bag. The impression is now safe and ready for further processing.



Switching to OPTIM

- Give all surfaces a good cleaning with soap and water before switching to OPTIM. After switching to OPTIM, surfaces need to be periodically rinsed with water (recommended no less than once per week)! Common practice for all chemical disinfectants is to periodically rinse and dry surfaces to maintain the surface integrity.
- No mixing of chemistries.