

Optim 1

SECTION 1. IDENTIFICATION

Product Identifier Optim 1
Recommended Use Ready to Use Disinfectant Cleaner.
Supplier SciCan Ltd, 1440 Don Mills road, Toronto, ON, M3B 3P9, 1-800-667-7733
Emergency Phone No. SciCan Ltd, 1-800-667-7733
SDS No. 053525

SECTION 2. HAZARD IDENTIFICATION

Classified according to Canada's Hazardous Products Regulations (WHMIS 2015)

GHS Label Elements

Signal Word: None
Hazard Pictogram: None
Hazard Statement(s): None
Precautionary Statement(s): None
Other Hazards

The product contains no substances which at their given concentration are considered to be hazardous to health.

SECTION 3. COMPOSITION/INFORMATION ON INGREDIENTS

Mixture:

Chemical Name	CAS No.	%	Other Identifiers
Hydrogen peroxide	7722-84-1	0.5	

SECTION 4. FIRST-AID MEASURES**First-aid Measures****Inhalation**

No specific first aid measures are required.

Skin Contact

No specific first aid measures are required.

Eye Contact

Flush with cool water. Remove contact lenses, if applicable, and continue washing. Obtain medical if irritation develops or persists.

Ingestion

No specific first aid measures are required.

SECTION 5. FIRE-FIGHTING MEASURES**Extinguishing Media****Suitable Extinguishing Media**

Not combustible. Use extinguishing agent suitable for surrounding fire.

Unsuitable Extinguishing Media

None known.

Specific Hazards Arising from the Product

None known.

Special Protective Equipment and Precautions for Fire-fighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

SECTION 6. ACCIDENTAL RELEASE MEASURES

Personal Precautions, Protective Equipment, and Emergency Procedures

Use the Personal Protective Equipment (PPE) recommended in Section 8.

Environmental Precautions

Before attempting clean-up, refer to hazard data. Prevent large spills from entering surface waterways.

Methods and Materials for Containment and Cleaning Up

Use water rinse for final cleanup.

SECTION 7. HANDLING AND STORAGE

Precautions for Safe Handling

Use good industrial hygiene practices in handling this material. FOR COMMERCIAL AND INDUSTRIAL USE ONLY.

Conditions for Safe Storage

Store in an area that is out of direct sunlight. Avoid storage at elevated temperatures. KEEP OUT OF REACH OF CHILDREN.

SECTION 8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Control Parameters

Chemical Name	ACGIH TLV®		OSHA PEL		AIHA WEEL	
	TWA	STEL	TWA	Ceiling	8-hr TWA	TWA
Hydrogen peroxide	1 ppm		1 ppm			

Appropriate Engineering Controls

No specific ventilation requirements.

Individual Protection Measures

Eye/Face Protection

Not required if product is used as directed.

Skin Protection

Not required if product is used as directed.

Respiratory Protection

Not required if product is used as directed.

SECTION 9. PHYSICAL AND CHEMICAL PROPERTIES

Basic Physical and Chemical Properties

Appearance	Clear colourless liquid.
Odour	Almond
Odour Threshold	Not available
pH	2.0 - 2.5
Melting Point/Freezing Point	Not available (melting); Not available (freezing)
Initial Boiling Point/Range	Not available
Flash Point	Not available
Evaporation Rate	Not available
Flammability (solid, gas)	Not applicable
Upper/Lower Flammability or Explosive Limit	Not available (upper); Not available (lower)
Vapour Pressure	Not available
Vapour Density (air = 1)	Not available
Relative Density (water = 1)	1.006 at 20 °C
Partition Coefficient, n-Octanol/Water (Log Kow)	Not available
Auto-ignition Temperature	Not available
Viscosity	1.09 centistokes at 20 °C

SECTION 10. STABILITY AND REACTIVITY

Reactivity

Not reactive.

Chemical Stability

This product is stable.

Possibility of Hazardous Reactions

None known.

Conditions to Avoid

Prolonged exposure to high temperatures.

Incompatible Materials

Do not mix with concentrated bleach products.

Hazardous Decomposition Products

None known.

SECTION 11. TOXICOLOGICAL INFORMATION

Likely Routes of Exposure

Inhalation; skin contact; eye contact; ingestion.

Acute Toxicity

LC50 (Inhalation): >2.08 mg/L

LD50 (Oral): >5000 mg/kg

LD50 (Dermal): >5050 mg/kg

Skin Corrosion/Irritation

Not classified under GHS criteria.

Serious Eye Damage/Irritation

Not classified under GHS criteria.

STOT (Specific Target Organ Toxicity) - Single Exposure

Inhalation

Not classified under GHS criteria.

Skin Absorption

Not classified under GHS criteria.

Ingestion

Not classified under GHS criteria.

STOT (Specific Target Organ Toxicity) - Repeated Exposure

Not classified under GHS criteria.

Respiratory and/or Skin Sensitization

Not classified under GHS criteria.

Carcinogenicity

Not classified under GHS criteria.

Reproductive Toxicity

Development of Offspring

Not classified under GHS criteria.

Sexual Function and Fertility

Not classified under GHS criteria.

Effects on or via Lactation

Not classified under GHS criteria.

Germ Cell Mutagenicity

Not classified under GHS criteria.

Interactive Effects

None known.

SECTION 12. ECOLOGICAL INFORMATION

This section is not required by WHMIS.

SECTION 13. DISPOSAL CONSIDERATIONS

Disposal Methods

Review federal, provincial and local government requirements prior to disposal.

SECTION 14. TRANSPORT INFORMATION

Not regulated under Canadian TDG regulations. Not regulated under US DOT Regulations.

Special Precautions Not applicable

Other Information Not regulated under IMO/IMDG

SECTION 15. REGULATORY INFORMATION

Safety, Health and Environmental Regulations

Canada

Domestic Substances List (DSL) / Non-Domestic Substances List (NDSL)

All ingredients are listed on the DSL/NDSL.

Additional Canadian Regulatory Lists

WHMIS exempt - DIN

SECTION 16. OTHER INFORMATION

HMIS Rating Health - 0 Flammability - 0 Physical Hazard - 0

SDS Prepared By Virox Technologies Inc.

Phone No. (800) 387-7578

Date of Preparation February 18, 2020

Additional Information

For an updated SDS please contact the supplier/manufacturer listed on the first page of this document. Information contained herein was obtained from sources considered technically accurate and reliable. While every effort has been made to ensure full disclosure of product hazards, in some cases data is not available and is so stated. Since conditions of actual product use are beyond control of the supplier, it is assumed that users of this material have been fully trained according to the requirement of all applicable legislation and regulatory instruments. No warranty, expressed or implied, is made and manufacturer/supplier will not be liable for any losses, injuries or consequential damages which may result from the use of or reliance on any information contained in this document. The contents of this document have been prepared in accordance with the Canadian Hazardous Products Regulations (WHMIS 2015).

Product Identifier: Optim 1

Date of Preparation: February 18, 2020

SDS No.: 053525

Page 04 of 04